

OMNI Tips™ plastic disposable generator probes

The OMNI Tips' patented design eliminates the hassle of disassembly and maintenance of traditional stainless steel generators.

Proprietary OMNI Tips™ Probes are available in both soft tissue (deaggregation) and hard tissue (frozen) versions. Engineered from durable plastics, OMNI Tips™ are reusable (autoclavable) and are so economical they can be disposed of after each use to eliminate cross contamination.

The outer tube of OMNI Tips™ is made of polycarbonate. The inner shafts are made of Ultem 1000, a polyetherimide. The OMNI Tips' patented design eliminates the hassle of disassembly and maintenance of traditional stainless steel generators. Treatable volume - 250 µL to 30 mL.



Features

- **Economical:** <1% of the cost of stainless steel probes
- **Eliminates Cross Contamination**
- **Ideal for PCR, RNA, and DNA sample preparation**
- **Maintenance Free**
- **Sterile OMNI Tips™ are Electron Beam treated to ensure they are free of contaminants.**

"With OMNI Tips™ I can be certain that all the tissue sample has been homogenized; there is no sample loss. I found that the OMNI Tips™ hard tissue version worked as well as the stainless steel probe."

— Rob T., National Institutes of Health

NEW Sterile OMNI Tips™ now available

A selection of our best-selling plastic probes:



SOFT TISSUE PLASTIC GENERATOR PROBES

	ORDER NO.
OMNI Tips™ Plastic Generator Probes, Soft Tissue 25 pack	30750
OMNI Tips™ Plastic Generator Probes, Soft Tissue 50 pack	32750
OMNI Tips™ Plastic Generator Probes, Soft Tissue 100 pack	34750

HARD TISSUE PLASTIC GENERATOR PROBES

	ORDER NO.
OMNI Tips™ Plastic Generator Probes, Hard Tissue 25 pack	30750H
OMNI Tips™ Plastic Generator Probes, Hard Tissue 50 pack	32750H
OMNI Tips™ Plastic Generator Probes, Hard Tissue 100 pack	34750H
Sterile OMNI Tips™ Generator Probes, Hard Tissue 25 pack	30750H-S

Note: A1000SB adapter is required for OMNI Tip™ use on OMNI GLH, GLH 850, PDH and TH (with black nose ring).

Note: 12-120 OMNI Tip™ adapter is required for OMNI Tip™ use on the OMNI TH115 and TH220 (with blue nose ring).



935 Cobb Place Blvd
Kennesaw, GA 30144
Phone: 770.421.0058 • Fax: 770.421.0206
www.omni-inc.com

Rev. F